

**NovaPED**

by **schein**

**PU foot orthotics**




## NOVAPED PU FOOT ORTHOTICS – BASIS

Novaped PU foot orthotics are made from a high-quality PU foam component. The structure of the foam is flexible, has a cushioning effect and relieves pressure. The resetting properties of the material ensure a reliable effect. The foot orthotics are of approx. 35° Shore A and in their basic version have an internal stiffening core. The orthotics are available with or without a meta pad.

Thanks to the raw material, the Novaped PU can be easily processed by the healthcare professional.

## Base module PU



### NOVAPED PU FOOT ORTHOTICS

- optional with or without meta pad
- internal stiffening core
- Basic module PU approx. 35° Shore A

Article no.	description	size
670100000	Foot orthotic with meta pad	35 - 50
670100300	Foot orthotic without meta pad	35 - 50



## NOVAPED PU FOOT ORTHOTICS – WITH TOP COVER

In combination with the new OnDry plus cover material, the wearing comfort of the basic module PU is further increased:

- Moisture regulation
- hypoallergenic functions
- roughened, velvety feel

The variant with 1.5 mm multiform orange offers additional cushioning comfort.


These two items can be provided with your logo on request.

### Base module PU with top cover

#### NOVAPED PU FOOT ORTHOTICS

- with meta pad, internal stiffening core
- Base module PU approx. 35° Shore A
- Top cover OnDry plus jeans

Article no.	description	size
670110347	Foot orthotic with meta pad, top cover OnDry plus jeans	35 - 50




### Base module PU with top layer and top cover

#### NOVAPED PU FOOT ORTHOTICS

- with meta pad, internal stiffening core
- Base module PU approx. 35° Shore A
- Multiform orange 1.5 mm approx. 30° Shore A
- Top cover OnDry plus melange

Article no.	description	size
670160363	Foot orthotic with meta pad, Multiform orange 1.5 mm, top cover OnDry plus melange	35 - 50



## NOVAPED PU FOOT ORTHOTICS – 2-LAYERS

This model is based on the PU foot orthotic and has additionally a grey Multiform cover layer in 1.5 mm thickness with a hardness of approx. 30° Shore A.

You can choose between different cover fabrics and optionally have your logo printed on the heel.

## Base module PU individual



### NOVAPED PU FOOT ORTHOTICS

- optional with or without meta pad
- Internal stiffening core
- Base module PU approx. 35° Shore A
- Multiform grey 1.5 mm approx. 30° Shore A

Article no.	description	size
670130150	Foot orthotic with meta pad, Multiform grey 1.5 mm	35 - 50
670130300	Foot orthotic without meta pad, Multiform grey 1.5 mm	35 - 50



## NOVAPED PU VARIO MODUL

The Novaped Vario Modul foot orthotic is characterized by different features.

Preset zones of loading and unloading:

- Cut-out of the base five
- Cushioning of the first and fifth metatarsal head

Individual surface design of the foot orthotic by the health care professional.

## Novaped PU Vario Modul



### NOVAPED PU VARIO MODUL

- Base module PU approx. 50° Shore A, top layer PU approx. 30° Shore A
- adaptable longitudinal arch wedge PU approx. 30° Shore A
- „PiGreco by schein“ microfibre top cover, blue, loose enclosed



Article no.	description	size
675110	Foot orthotic wide	39/40 - 47/48
675210	Foot orthotic narrow	35/36 - 41/42



03.0997.02.07.20



099993213

Schein Orthopädie Service KG  
Hildegardstr. 5  
42897 Remscheid  
Germany  
Tel. +49 2191 910-200  
Fax +49 2191 910-100  
remscheid@schein.de  
www.schein.de



— Germany since 1879 —